

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution
Product code(s): 04-182
Batch Number(s): 32181029
Date Manufactured: 15 OCT 18
Expiry Date: 24 SEP 19
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.05
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	11780
Osmolality(mOsm/kg)	N/A	N/A

NHS BT Spec
Manufactured to SPN870/2 Dated 18/08/16

Preservatives used: Chloramphenicol and Neomycin sulphate

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2-8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 13 NOV 18

Signed for Certificate Check: 

Date: 13 Nov 18