

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsevers Solution  
**Product code(s):** 04-182  
**Batch Number(s):** 32181030  
**Date Manufactured:** 29 OCT 18  
**Expiry Date:** 24 SEP 19  
**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	6.98
Conductivity( $\mu\text{s}/\text{cm}$ )	11500 - 13000	11980
Osmolality(mOsm/kg)	N/A	N/A

NHSBT Spec  
SPN870/2, Dated 18/08/16

**Preservatives used:** Chloramphenicol and Neomycin sulphate

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2 - 8°C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 13 NOV 18

**Signed for Certificate Check:** 

**Date:** 13 Nov 18