

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Phosphate Buffered Saline Concentrate (PBSC)  
**Product code(s):** 04-007, 04-008 and 04-009  
**Batch Number(s):** 04180710  
**Date Manufactured:** 04 Jul 18  
**Expiry Date:** 01 Oct 19  
**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	6.95-7.05	6.95
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	285-295	289

**Preservatives used:** Sulfamethoxazole, Trimethoprim, Sodium azide (0.01% w/v)

**Sterility:** N/A

**Storage:** Ambient

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

15 Nov 18

**Signed for Certificate Check:**



**Date:**

15 Nov 18