

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsevers Solution

**Product code(s):** 05-182

**Batch Number(s):** 99181001

**Date Manufactured:** 08 Oct 18

**Expiry Date:** 25 Sep 19

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.80-7.30	6.97
Conductivity( $\mu\text{s}/\text{cm}$ )	11500-13000	11500
Osmolality(mOsm/kg)	335-360	352

**Preservatives used:** Chloramphenicol, Neomycin sulphate and Levofloxacin

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2 - 8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

19 Nov 18.

**Signed for Certificate Check:**



**Date:**

19 Nov 18.