

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution
Product code(s): 05-182
Batch Number(s): 99181102
Date Manufactured: 12 Nov 18
Expiry Date: 12 Jan 20
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80-7.30	7.12
Conductivity($\mu\text{s}/\text{cm}$)	11500-13000	11930
Osmolality(mOsm/kg)	335-360	351

Preservatives used: Chloramphenicol, Neomycin sulphate and Levofloxacin

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date:

20 Nov 18

Signed for Certificate Check:



Date:

20 Nov 18