

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Frozen Cell Recovery Solution (FCRS)

Product code(s): 04-161 Bag

Batch Number(s): 41180918

Date Manufactured: 19 Sep 18

Expiry Date: 04 Apr 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.40 - 6.80	6.59
Conductivity(μs/cm)	44000 - 52000	48136
Osmolality(mOsm/kg)	1080-1120	1109

NHS BT Spec

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: EDTA (anti-fungal)

Sterility: Product is filtered through a 0.2 μ m filter into gamma irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 11 Oct 18

Signed for Certificate Check:



Date: 11 OCT 18.