

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cryoprotectant Reagent (CPR)

Product code(s): 04-213

Batch Number(s): 61181112

Date Manufactured: 21 Nov 18

Expiry Date: 31 Mar 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	4.00-6.00	5.57
Conductivity($\mu\text{s}/\text{cm}$)	4200-4850	4425
Osmolality(mOsm/kg)	880-970	949

Preservatives used: N/A

Sterility: Terminally filtered with 0.2 μm filter in a clean room and tested for sterility

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 03 Dec 18

Signed for Certificate Check:



Date: 03 DEC 18