

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: LISS (Low Ionic Strength Saline)

Product code(s): 04-016

Batch Number(s): 03181128

Date Manufactured: 30 Nov 18

Expiry Date: 29 Aug 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.50-6.90	6.57
Conductivity($\mu\text{s}/\text{cm}$)	3500-3800	3656
Osmolality(mOsm/kg)	285-295	289

Preservatives used: Sodium azide (0.01% [w/v])

Sterility: 0.2 μ Filter sterilised in a clean room

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 03 Dec 18

Signed for Certificate Check:



Date: 03 Dec 18