

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Decellularised Dermis Cryopreservation Medium
Product code(s): 04-351
Batch Number(s): 83181113
Date Manufactured: 27 Nov 18
Expiry Date: 30 Sep 20
Appearance: Clear Liquid
Quality Control:

	Specification	This batch result
pH	7.10-7.50	7.28
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

NHS BT Spec
Meets requirements set of in SPN430/3

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter in a clean room and tested for sterility.

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 03 Dec 18

Signed for Certificate Check:



Date: 03 Dec 18