

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** PBSS (PHOSPHATE BUFFERED SEROLOGY SALINE )  
**Product code(s):** 04-004  
**Batch Number(s):** 01181288  
**Date Manufactured:** 05 Dec 18  
**Expiry Date:** 04 Dec 20  
**Appearance:** Clear Liquid  
**Quality Control:**

|                                         | Specification | This batch result |
|-----------------------------------------|---------------|-------------------|
| pH                                      | 6.95-7.05     | 6.98              |
| Conductivity( $\mu\text{s}/\text{cm}$ ) | N/A           | N/A               |
| Osmolality(mOsm/kg)                     | 285-295       | 292               |

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** Sodium azide (0.01% [w/v])

**Sterility:** Product is filled into gamma irradiated bags

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 06 Dec 18

Signed for Certificate Check: 

Date: 06 Dec 18