

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Fixol (Fixative reagent for Fetal Cell Detection Kit)
Product code(s): 05-051
Batch Number(s): 51181146
Date Manufactured: 06 Nov 18
Expiry Date: 16 Oct 20
Appearance: Colourless transparent liquid

Quality Control:

	Specification	This batch result
pH	8.50 - 9.50	8.70
Conductivity($\mu\text{s}/\text{cm}$)	1750 - 2150	1912
Osmolality(mOsm/kg)	N/A	N/A

The performance of the above components has been assessed by Kleihauer testing of negative and positive control samples. All tests met the performance specifications for this product and passed the assessment. [Demonstration of fetal hemoglobin in erythrocytes of a blood smear]. KLEIHAUER E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

Preservatives used: Sodium azide (0.024% w/v), Alcohol based

Sterility: N/A
Storage: Room Temperature
Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 11 Dec 18

Signed for Certificate Check:



Date: 11 Dec 18

CERTIFICATE OF ANALYSIS