

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Low Ionic Strength Saline (LISS)

Product code(s): 04-016

Batch Number(s): 03181229

Date Manufactured: 13 Dec 18

Expiry Date: 29 Aug 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.50-6.90	6.77
Conductivity($\mu\text{s}/\text{cm}$)	3500-3800	3628
Osmolality(mOsm/kg)	285-295	285

Preservatives used: Sodium Azide (0.09% [w/v])

Sterility: N/A

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Wear appropriate PPE when handling this reagent.


Please download the safety data sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 17 DEC 18

Signed for Certificate Check:



Date: 17 Dec 18