

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Reagent Grade Water

**Product code(s):** 04-035

**Batch Number(s):** 18181298

**Date Manufactured:** 13 Dec 18

**Expiry Date:** 12 Dec 19

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	N/A	N/A
Conductivity( $\mu$ s/cm)	0 - 10	2.736
Osmolality(mOsm/kg)	N/A	N/A

**Preservatives used:** No Preservatives used

**Sterility:** Water is UV treated and processed into gamma irradiated bags, Although exposure to the environment is limited, we do not guarantee the sterility of the water.

**Storage:** Store at Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

There are no H&S implications from this product ;

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 17 Dec 18

**Signed for Certificate Check:** 

**Date:** 17 Dec 18