

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: PBSS AF (PHOSPHATE BUFFERED SEROLOGY SALINE AZIDE -FREE)

Product code(s): 04-103

Batch Number(s): 01171295

Date Manufactured: 20 Dec 18

Expiry Date: 19 Dec 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.95-7.05	6.99
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295	285

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: N/A

Sterility: Product is filled into gamma irradiated bags

Storage: Store at Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 03 JAN 19

Signed for Certificate Check:



Date: 03 Jan 19