

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Low Ionic Strength Saline (LISS)

**Product code(s):** 04-015, 04-016, 04-017

**Batch Number(s):** 03190130

**Date Manufactured:** 22 Jan 19

**Expiry Date:** 04 Apr 20

**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	6.50-6.90	6.84
Conductivity( $\mu\text{s}/\text{cm}$ )	3500-3800	3670
Osmolality(mOsm/kg)	285-295	292

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** Sodium azide (0.01% [w/v])

**Sterility:** N/A

**Storage:** Room temperature

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 25 Jan 19

**Signed for Certificate Check:**



**Date:** 25 JAN 19