

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** PBSS ( PHOSPHATE BUFFERED SEROLOGY SALINE )  
**Product code(s):** 04-004 - 04-005  
**Batch Number(s):** 01190345  
**Date Manufactured:** 15 Mar 19  
**Expiry Date:** 14 Mar 21  
**Appearance:** Clear Liquid  
**Quality Control:**

	Specification	This batch result
pH	6.95-7.05	6.95
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	285-295	292

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** Sodium azide (0.01% [w/v])

**Sterility:** Product is filled into gamma and Pinto irradiated bags

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

19 MAR 19

**Signed for Certificate Check:**



**Date:**

17 MAR 19