

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Skin Cryopreservation Medium  
**Product code(s):** 04-336  
**Batch Number(s):** 73190407  
**Date Manufactured:** 08 Apr 19  
**Expiry Date:** 10 Jun 20  
**Appearance:** Clear Red Solution

**Quality Control:**

	Specification	This batch result
pH	7.10-7.50	7.22
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN311/2, effective 29/03/17.

**Preservatives used:** N/A

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter and is tested for Micro Sterility

**Storage:** 2-25° C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 15 Apr 19

**Signed for Certificate Check:**



**Date:** 15 Apr 19