

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Skin Cryopreservation Medium

Product code(s): 04 - 336

Batch Number(s): 73190408

Date Manufactured: 25 Apr 19

Expiry Date: 01 Sep 20

Appearance: Clear Red Solution

Quality Control:

	Specification	This batch result
pH	7.10 - 7.50	7.30
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN311/2, effective 29/03/17.

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter and is tested for Micro Sterility

Storage: 2-25° C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 30 APR 19

Signed for Certificate Check: 

Date: 30 Apr 19