

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Polyagglutination Kit
Product code(s): 03-061
Batch Number(s): 79190514
Date Manufactured: 01 May 19
Expiry Date: 16 Apr 21
Appearance: Clear Yellow/Amber Liquid(s),
Kit made up of 4 lectin extracts: 1. Arachis hypogea - 14190409, 2. Glycine soja - 15190308
3. Salvia Sclarea - 79190403, 4. Dolichos biflorus - 27190312

Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | N/A | N/A |
| Conductivity($\mu\text{s}/\text{cm}$) | N/A | N/A |
| Osmolality(mOsm/kg) | N/A | N/A |

No specific guidelines are laid down for a polyagglutination kit, but the kit is expected to show a pattern of agglutination with affected human red cells which will classify them as being T, Tn, Tk or Cad. It may be possible to further determine Th and Tx affected red cells.

Preservatives used: Sodium azide (0.01% [w/v])

Sterility: N/A

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 02 May 19

Signed for Certificate Check:



Date: 02 May 19