

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Liss 2 (Liss Additive Reagent)

Product code(s): 04 - 070

Batch Number(s): 34190517

Date Manufactured: 08 May 19

Expiry Date: 04 Jun 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.70 - 7.10	7.08
Conductivity($\mu\text{s}/\text{cm}$)	3000 - 3600	3260
Osmolality(mOsm/kg)	290 - 310	305

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: Sodium azide (0.07% [w/v])

Sterility: Product is filtered through a 0.2 μm filter in a clean room and tested for sterility

Storage: 2-8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 14 MAY 19

Signed for Certificate Check: 

Date: 14 MAY 19