

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Fixol (Fixative reagent for Fetal Cell Detection Kit)  
**Product code(s):** 05-051  
**Batch Number(s):** 51190453  
**Date Manufactured:** 16 Apr 19  
**Expiry Date:** 24 Mar 21  
**Appearance:** Clear transparent liquid  
**Quality Control:**

	Specification	This batch result
pH	8.50 - 9.50	9.19
Conductivity( $\mu\text{s}/\text{cm}$ )	1750 - 2150	1944
Osmolality(mOsm/kg)	N/A	N/A

This product has been assessed and Passed using the Kleihaur testing of negative and positive control samples. (Demonstration of fetal hemoglobin in erythrocytes of a blood smear). KLEIHAUR E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

**Preservatives used:** Methanol based product

**Sterility:** N/A

**Storage:** Ambient (Room Temperature)

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 07 May 19

**Signed for Certificate Check:** 

**Date:** 07 May 19