

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Anti-H (Ulex Europaeus Seeds)

Product code(s): 03-007

Batch Number(s): 10190502

Date Manufactured: 15 May 19

Expiry Date: 14 May 21

Appearance: Clear Yellow Liquid

Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	290-1000	324

Preservatives used: Sodium azide (0.01% [w/v])

Sterility: N/A

Storage: 2-8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 23 May 19

Signed for Certificate Check:



Date: 23 MAY 19