

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cambridge Antibiotic Solution

Product code(s): 04-301

Batch Number(s): 70191101

Date Manufactured: 26 Nov 19

Expiry Date: 30 Sep 20

Appearance: Red Suspension

Quality Control:

	Specification	This batch result
pH	7.10 - 7.50	7.19
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	280 - 300	292

Manufactured to SPN 310/3 effective 07 May 19

Preservatives used: Gentamycin Sulphate, Imipenem, Polymixin B, Vancomycin and Nystatin

Sterility: Product is filtered through a 0.2 μm filter prior to addition of Nystatin, in a clean room and filled into gamma irradiated bottles

Storage: <-40°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 02 Dec 19

Signed for Certificate Check: 

Date: 02 Dec 19