

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Amniotic Membrane Cryopreservation Medium
Product code(s): 04-326
Batch Number(s): 72200101
Date Manufactured: 15 Jan 20
Expiry Date: 01 Jan 21
Appearance: Clear Red Solution

Quality Control:

	Specification	This batch result
pH	7.10 - 7.50	7.50
Conductivity(μ s/cm)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN431/2 Effective 06 Jul 17

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μ m filter into gamma irradiated bottles


Storage: 0 - 30 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 20 Jan 2020

Signed for Certificate Check: 

Date: 20 JAN 2020