

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Detergent Buffer Base

**Product code(s):** 04-334

**Batch Number(s):** 77200201

**Date Manufactured:** 19 Feb 2020

**Expiry Date:** 31 Oct 2020

**Appearance:** Clear Solution

**Quality Control:**

	Specification	This batch result
pH	7.8 - 8.2	8.08
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN 321/2 Effective 05Apr17

**Preservatives used:** N/A

**Sterility:** Product is ran through a 0.2 $\mu\text{m}$  filter into gamma irradiated bottles

**Storage:** Room Temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 24 Feb 2020

**Signed for Certificate Check:**



**Date:** 24 Feb 2020