

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Hypotonic Buffer Base
Product code(s): 04 - 333
Batch Number(s): 76200201
Date Manufactured: 19 Feb 2020
Expiry Date: 28 Jan 2021
Appearance: Clear Liquid
Quality Control:

	Specification	This batch result
pH	7.80 - 8.20	8.10
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN320/2 Effective 29 Mar 17

Preservatives used: N/A

Sterility: Product is ran through a 0.2 μm filter and filled into Gamma Irradiated Bottles


Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 25 Feb 2020

Signed for Certificate Check: 

Date: 25 FEB 2020