

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Soutlion

Product code(s): 05 - 182

Batch Number(s): 99200201

Date Manufactured: 25 Feb 2020

Expiry Date: 24 Feb 2021

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	6.97
Conductivity(μ s/cm)	11500 - 13000	11753
Osmolality(mOsm/kg)	335 - 360	350

Preservatives used: Chloramphenicol, Neomycin Sulphate and Levofloxacin

Sterility: Product is filtered through a 0.2 μ m filter into Gamma Irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 02 MAR 2020

Signed for Certificate Check: 

Date: 02 MAR 2020