

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Low Ionic Strength Saline (LISS)

**Product code(s):** 04-016, 04-017 and 04-123

**Batch Number(s):** 03200301

**Date Manufactured:** 09 MAR 2021

**Expiry Date:** 08 MAR 2021

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.50 - 6.90	6.86
Conductivity( $\mu\text{s}/\text{cm}$ )	3500 - 3800	3558
Osmolality(mOsm/kg)	285 - 295	290

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** Sodium azide (0.1% [w/v])

**Sterility:** Product is filled into gamma irradiated bags and bottles

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 13 MAR 2020

**Signed for Certificate Check:**



**Date:** 13 MAR 2020