

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Tris Buffered Isotonic Saline (TBIS)

Product code(s): 04 - 142

Batch Number(s): 37200301

Date Manufactured: 25 Mar 2020

Expiry Date: 24 Mar 2021

Appearance: Clear Fluid

Quality Control:

	Specification	This batch result
pH	7.3 - 7.6	7.56
Conductivity($\mu\text{s}/\text{cm}$)	14000 - 16500	14250
Osmolality(mOsm/kg)	285 - 300	298

Manufactured to SPN 437/1 Effective 11May12

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into Gamma Irradiated Bottles.

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *S.J. Taylor*

Date: 30MAR2020

Signed for Certificate Check: *T Waller*

Date: 30MAR 2020