

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cardiovascular and Meniscus Cryopreservation Medium
Product code(s): 04-311
Batch Number(s): 71200401
Date Manufactured: 07 APR 2020
Expiry Date: 31 Dec 2020
Appearance: Clear Red Solution

Quality Control:

	Specification	This batch result
pH	7.0 - 7.5	7.39
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN323/3 Effective 24 Sep 18

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bottles

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 09 APR 2020

Signed for Certificate Check:



Date: 14 Apr 2020