

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cryoprotectant Reagent (CPR)

Product code(s): 04 - 213

Batch Number(s): 61200401

Date Manufactured: 07 APR 2020

Expiry Date: 09 APR 2021

Appearance: Clear, Viscous Liquid

Quality Control:

	Specification	This batch result
pH	5.80 - 6.70	6.41
Conductivity($\mu\text{s}/\text{cm}$)	4200 - 4850	4260
Osmolality(mOsm/kg)	880 - 970	931

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into Gamma Irradiated Bags

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:

S. J. Taylor

Date: 20 APR 2020

Signed for Certificate Check:

Waller

Date: 20 APR 2020