

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution

Product code(s): 04 - 182, 04 - 182HE

Batch Number(s): 32200601

Date Manufactured: 02 JUNE 2020

Expiry Date: 01 JUNE 2021

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.16
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	12443
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN 870/2.1 Dated 06.11.18

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 09 June 2020

Signed for Certificate Check:



Date: 09 JUN 2020