

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution

Product code(s): 05 - 182

Batch Number(s): 99200601

Date Manufactured: 24 Jun 2020

Expiry Date: 22 Jun 2021

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.06
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	11657
Osmolality(mOsm/kg)	335 - 360	359

Preservatives used: Chloramphenicol, Neomycin Sulphate and Levofloxacin

Sterility: Product is filtered through 0.2 μm filter into Gamma Irradiated Bottles

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *S.J. Taylor*

Date: 02 Jul 2020

Signed for Certificate Check: *T. Waller*

Date: 02 Jul 2020