

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution
Product code(s): 05 - 182
Batch Number(s): 99200602
Date Manufactured: 25 Jun 2020
Expiry Date: 22 Jun 2021
Appearance: Clear Liquid

Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | 6.80 - 7.30 | 7.14 |
| Conductivity($\mu\text{s}/\text{cm}$) | 11500 - 13000 | 11560 |
| Osmolality(mOsm/kg) | 335 - 360 | 351 |

Preservatives used: Chloramphenicol, Neomycin Sulphate and Levofloxacin

Sterility: Product is filtered through a 0.2 μm filter into Gamma Irradiated Bags

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 02 Jul 2020

Signed for Certificate Check: 

Date: 02 July 2020