

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Frozen Cell Recovery Solution (FCRS)
Product code(s): 04 - 161 BAG
Batch Number(s): 41200601
Date Manufactured: 30 Jun 2020
Expiry Date: 08 Sep 2021
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.40 - 6.80	6.54
Conductivity(μ s/cm)	44000 - 52000	46980
Osmolality(mOsm/kg)	1080 - 1120	1105

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μ m filter into Gamma Irradiated Bags

Storage: 2 - 8°C

Country of Origin: United Kingdom


Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 06 Jul 2020

Signed for Certificate Check: 

Date: 06 July 2020