

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Isotonic Saline  
**Product code(s):** 04 - 025  
**Batch Number(s):** 66200701  
**Date Manufactured:** 13 July 2020  
**Expiry Date:** 12 May 2021  
**Appearance:** Clear, Colourless Liquid  
**Quality Control:**

	Specification	This batch result
pH	N/A	N/A
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	285 - 295	292

**Preservatives used:** N/A

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into Gamma Irradiated Bottles


**Storage:** Room Temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 21 Jul 2020

**Signed for Certificate Check:** 

**Date:** 21 July 2020