

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Phosphate Buffered Serology Saline Concentrate  
**Product code(s):** 04-007, 04-008  
**Batch Number(s):** 04200701  
**Date Manufactured:** 21 Jul 20  
**Expiry Date:** 12 May 21  
**Appearance:** Clear Liquid  
**Quality Control:**

	Specification	This batch result
pH	6.95 - 7.05	6.95
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	285 - 295	289

**Preservatives used:** Sulfamethoxazole, Trimethoprim, Sodium Azide (0.01% [w/v])

**Sterility:** N/A  
**Storage:** Room temperature  
**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 27 July 2020

Signed for Certificate Check: 

Date: 27 JUL 2020