

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cambridge Antibiotic Solution

Product code(s): 04-301

Batch Number(s): 70200701

Date Manufactured: 16JUL2020

Expiry Date: 16APR21

Appearance: Red Suspension

Quality Control:

| | Specification | This batch result |
|---------------------|---------------|-------------------|
| pH | 7.1 - 7.5 | 7.22 |
| Osmolality(mOsm/kg) | 280 - 300 | 299 |

Manufactured to SPN 310/3 effective 07 May 19

Preservatives used: Gentamycin Sulphate, Imipenem, Polymixin B, Vancomycin, Nystatin

Sterility: Product is filtered through a 0.2µm filter prior to addition of Nystatin in a clean room and filled into Gamma Irradiated bottles.

Storage: <-40 °C


Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 31 Jul 2020

Signed for Certificate Check:



Date: 31 July 2020