

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsevers Solution

**Product code(s):** 04-182

**Batch Number(s):** 32200801

**Date Manufactured:** 27 AUG 2020

**Expiry Date:** 12 MAY 2021

**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.14
Conductivity( $\mu\text{s}/\text{cm}$ )	11500 - 13000	12400

Manufactured to SPN 870/2.1 Dated 18.08.16

**Preservatives used:** Chloramphenicol and Neomycin Sulphate

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2 - 8°C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes


Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 01 SEP 2020

**Signed for Certificate Check:**



**Date:** 01 Sep 2020