

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Fixol (Fixative Reagent for Fetal Cell Detection Kit)
Product code(s): 05 - 051
Batch Number(s): 51200801
Date Manufactured: 07 Sep 2020
Expiry Date: 01 Sep 2021
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	8.50 - 9.50	8.62
Conductivity($\mu\text{s}/\text{cm}$)	17500 - 21500	18860
Osmolality(mOsm/kg)	N/A	N/A

The performance of the above components has been assessed by Kleihauer testing of negative and positive control samples. All tests met the performance specifications for this product and passed the assessment. [Demonstration of fetal hemoglobin in erythrocytes of a blood smear]. KLEIHAUER E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

Preservatives used: Sodium Azide (0.024% w/v)

Sterility: N/A

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 23 Sep 2020

Signed for Certificate Check:



Date: 23 Sep 2020

Author: Thomas Waller
Technical Reviewer: Manuara Mahmood
QA Reviewer: Jahira Jangaria

SER-DOC-32 v3.0
Authorised Date: 03 Dec 19
Next Review Date: 03 Dec 21