

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Fetal Cell Detection Kit (FCD Kit)
Product code(s): 05 - 050
Batch Number(s): 05200801 A/B
Date Manufactured: 20 Aug 2020
Expiry Date: 03 Apr 2021
Appearance: Eluting Reagent - Brown/Black Opaque Liquid
Counterstain - Pink Translucent Liquid

Quality Control:

	Specification	This batch result
pH	Eluting - 1.80 - 2.10	Eluting - 1.91
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

The performance of the above components has been assessed and met the standard set by Kleihauer testing of positive control samples. [Demonstration of fetal hemoglobin in erythrocytes of a blood smear].
KLEIHAUER E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

Preservatives used: Eluting - Alcohol Based, Counterstain - Sodium Azide (0.10% w/v)

Sterility: N/A

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 23 Sep 2020

Signed for Certificate Check:



Date: 23 Sep 2020