

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: AmnioCryo Medium (Amniotic Membrane Cryopreservation)
Product code(s): 04-326
Batch Number(s): 72200901
Date Manufactured: 22 Sep 2020
Expiry Date: 12 Jan 2022
Appearance: Clear Red Solution
Quality Control:

	Specification	This batch result
pH	7.10-7.50	7.43

Manufactured to SPN431/2 – Effective date 06 July 2017.

Preservatives used: N/a

Sterility: Product is filtered through a 0.2µm filter and filled into gamma irradiated bottles..

Storage: Ambient.


Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 02 Oct 2020

Signed for Certificate Check:

Date: 02 Oct 2020