

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Phosphate Buffered Seroloy Saline (Azide Free)
Product code(s): 04 - 047
Batch Number(s): 25201001
Date Manufactured: 07 Oct 2020
Expiry Date: 31 Dec 2021
Appearance: Clear Liquid
Quality Control:

	Specification	This batch result
pH	7.25 - 7.35	7.31
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285 - 295	289

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into Gamma Irradiated Bottles


Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 12 OCT 2020

Signed for Certificate Check: 

Date: 12 Oct 2020