

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cardiovascular and Meniscus Cryopreservation Medium
(Heart Valve Cryopreservation Medium))

Product code(s): 04 - 311

Batch Number(s): 71201001

Date Manufactured: 05 Oct 2020

Expiry Date: 31 Jan 2022

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.00 - 7.50	7.49
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN 323/3 (SER SPEC 10) Effective 24.09.18

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into Gamma Irradiated Bottles

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 12 Oct 2020

Signed for Certificate Check: 

Date: 12 Oct 2020