

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Reagent grade water

Product code(s): 04 - 035

Batch Number(s): 18201101

Date Manufactured: 19 Nov 2020

Expiry Date: 18 Nov 2022

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
Conductivity($\mu\text{s}/\text{cm}$)	< 10	4.475
Osmolality(mOsm/kg)	< 10	1

Preservatives used: N/a

Sterility: Product is dispensed into gamma irradiated bags

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 25 Nov 2020

Signed for Certificate Check: 

Date: 26 Nov 2020