

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cambridge Antibiotic Solution

Product code(s): 04-301

Batch Number(s): 70201001

Date Manufactured: 08 Oct 2020

Expiry Date: 31 Jul 2021

Appearance: Red Suspension

Quality Control:

	Specification	This batch result
pH	7.1 - 7.5	7.12
Osmolality(mOsm/kg)	280 - 300	297

Manufactured to SPN 310/3 effective 07 May 19

Preservatives used: Gentamycin Sulphate, Imipenem, Polymixin B, Vancomycin, Nystatin

Sterility: Product is filtered through a 0.2µm filter prior to addition of Nystatin in a clean room and filled into Gamma Irradiated bottles.

Storage: <-40 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 04 DEC 2020

Signed for Certificate Check:



Date: 04 Dec 2020