

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Isotonic Saline
Product code(s): 04-025
Batch Number(s): 66210101
Date Manufactured: 05 Jan 21
Expiry Date: 04 Jan 22
Appearance: Clear, Colourless Liquid
Quality Control:

	Specification	This batch result
Osmolality(mOsm/kg)	285 - 295	288

Preservatives used: N/A

Sterility: Product is filtered through a 0.2µm filter into Gamma Irradiated Bottles

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date:

13 JAN 21

Signed for Certificate Check:



Date:

13 JAN 21