

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: 0.5M Hydrochloric acid

Product code(s): 04-356

Batch Number(s): 84210101

Date Manufactured: 12 Jan 2021

Expiry Date: 13 Feb 2022

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
Osmolality(mOsm/kg)	920-1020	989

Manufactured to SPN 435/1 – Effective date 28 Mar 2012

Preservatives used: N/a

Sterility: Product is filtered through a 0.2µm filter into gamma irradiated bottles.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 15 JAN 21

Signed for Certificate Check: 

Date: 15 Jan 21