

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Red Cell Support Solution (RCSS)

**Product code(s):** 04-131Bag

**Batch Number(s):** 35210101

**Date Manufactured:** 18 Jan 21

**Expiry Date:** 31 Dec 21

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.80 - 7.20	7.07
Conductivity( $\mu\text{s}/\text{cm}$ )	3100 - 3300	3264
Osmolality(mOsm/kg)	300 - 315	306

**Preservatives used:** Antibiotics - Sulfamethoxazole and Trimethoprim

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2-8 $^{\circ}\text{C}$

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 25 JAN 21

**Signed for Certificate Check:**



**Date:** 25 Jan 21