

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Tris Buffered Isotonic Saline (TBIS)

Product code(s): 04 - 142 (Ser spec 16 / SPN 437.1)

Batch Number(s): 37210301

Date Manufactured: 04 Mar 2021

Expiry Date: 03 Mar 2022

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.3 – 7.6	7.39
Conductivity($\mu\text{s}/\text{cm}$)	14000 - 16500	14433
Osmolality(mOsm/kg)	285 – 300	299

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bottles

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 12 March 21

Signed for Certificate Check:



Date: 12 MAR 21